

UV400 UPGRADE OF THE UV PROTECTION



1.50 UV400

1.56 1.60 1.67 1.74 ALL INDEX

Transparent lens with
The same UV-protection
Function as the sunglasses



UV400

The protection
function
is not affected
by coating.

www.conantoptical.com

Conant[®]

The imported raw materials with UV protection function. It has long-term protection for UV light and the protection function is not affected by coating in monomer.

Conant UV400 lens makes use of advanced technology that blocks UV light, and the colorless raw material guarantee the transparency of the lens. At the same time, the lens does not use a coating to block the UV light, so the UV protection function is not affected by the coating.




Transparent lens with the same UV-protection function as the sunglasses

Sunglasses are one of the most popular UV blocking products. The UV400 lens is an all-weather lens with excellent UV protection function. It not only has the same UV blocking effect as standard sunglasses, but also can allow us to wear it in any environment.

1.50 UV400				
Process	Base Color	Process	Index	Abe Value
The monomer absorbs UV light	transparent	100% UV-protection	1.50	58

“ UV400 ”

“ Upgraded UV protection.
100% UV protection up to 400nm ”



At present, the UV blocking standard of most lenses is UV380, which can't reach the UV blocking standard of UV400. Conant's UV400 lens can achieve 100% UV protection for wave lengths up to 400nm. Conant's lens can protect your eyes from the damages of the harmful ultraviolet light.

